



REF: 225 001 60 Test
REF: 225 002 120 Test

Normal and abnormal commercial control serum of known concentrations should be analyzed with each run.

Performance Characteristics

Precision

Within run (Repeatability)

| | Total | | Direct | |
|--------------|---------|---------|---------|---------|
| | Level 1 | Level 2 | Level 1 | Level 2 |
| n | 20 | 20 | 20 | 20 |
| Mean (mg/dL) | 0.79 | 4.37 | 0.299 | 0.77 |
| SD | 0.016 | 0.18 | 0.016 | 0.057 |
| CV% | 2.03 | 4.12 | 5.35 | 7.4 |

Run to run (Reproducibility)

| | Total | | Direct | |
|--------------|---------|---------|---------|---------|
| | Level 1 | Level 2 | Level 1 | Level 2 |
| n | 20 | 20 | 20 | 20 |
| Mean (mg/dL) | 0.82 | 4.52 | 0.32 | 0.82 |
| SD | 0.02 | 0.17 | 0.023 | 0.062 |
| CV% | 2.44 | 3.76 | 7.19 | 7.56 |

Methods Comparison

A comparison between Spectrum Bilirubin and a commercial reagent of the same methodology was performed on 200 human sera. A correlation of 0.975 was obtained.

Sensitivity

When run as recommended, the sensitivity of this assay is 0.1 mg/dL (1.7 µmol/L) for total and 0.04 mg/dL (0.68 µmol/L) for direct bilirubin.

Linearity

The reaction is linear up to a total bilirubin concentration of 18 mg/dL (308 µmol/L) and a direct bilirubin concentration of 18 mg/dL (308 µmol/L). Specimens showing higher concentration should be diluted 1+4 with physiological saline and repeat the assay (result×5).

Interfering substances

Haemolysis

Avoid haemolysis since it interferes with the test.

Lipemia

Lipemic specimens interfere with the test.

Drugs

Theophylline and propranolol may cause artificially low total bilirubin levels.

Expected Values

Total Bilirubin

Adults and infants >1 month < 0.2-1.1 mg/dL (3.4-17 µmol/l)
Newborns premature (3-5 d) 10-14 mg/dL (171-239 µmol/l)

Newborns:

(3-5 d) 4.0 - 8.0 mg/dL (68-137 µmol/l)
(<48 h) 6.0 - 10.0 mg/dL (103-171 µmol/l)
(<24 h) 2.0-6.0 mg/dL (34-103 µmol/l)

Direct Bilirubin 0 - 0.3 mg/dL (0 - 5.1 µmol/L)

Spectrum Diagnostics does not interpret the results of a clinical laboratory procedure; interpretation of the results is considered the responsibility of qualified medical personnel. All indications of clinical significance are supported by literature references.

Analytical Range

Total bilirubin : 0.1 – 18 mg/dL (1.7 – 308 µmol/L)
Direct bilirubin : 0.04 – 18 mg/dL (0.68 – 308 µmol/L)

Waste Disposal

This product is made to be used in professional laboratories. Please consult local regulations for a correct waste disposal.

S56: dispose of this material and its container at hazardous or special waste collection point.

S57: use appropriate container to avoid environmental contamination.

S61: avoid release in environment. refer to special instructions/safety data sheets.

References

1. Balistreri WF, Shaw LM. Liver function. In: Tietz NW, ed. Fundamentals of clinical chemistry. 3rd ed. Philadelphia: WB Saunders; 1987:729-761.
2. Malloy HT, Evelyn KA. The determination of bilirubin with the photoelectric colorimetric method. J Biol Chem. 1937;119:481-490.
3. Tietz NW, ed. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders; 1995:268-273.



Egyptian Co for Biotechnology - Spectrum Diagnostics (S.A.E)

Obour city industrial area, block 20008 piece 19 A. Cairo. Egypt.

Tel: +202 4489 2248 - Fax: +202 4489 2247

www.spectrum-diagnostics.com

E-mail: info@spectrum-diagnostics.com



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



IFUFCC53

Rev.(6), 6/6/2021