

# SP-Multical

REF: 328 001 (1 vial ) 1 x 3 ml  
REF: 328 002 (5 vials) 5 x 3 ml

## Intended use

Multiparameter calibration serum for the standardization of quantitative tests on clinical chemical Analytes

## Description

SP-MULTICAL is an analyzed calibration material based on lyophilized human serum.

## Storage and stability

Unopened bottles have to be stored in refrigerator at 2 – 8 °C and are stable up to the expiry date printed on the labels.

After reconstitution the tightly closed calibrator can be used within the times and temperatures given below:

### General:

18 – 25°C => 8 hours  
2 – 8°C => 2 days  
< – 25°C => 30 days

### Bilirubin (in the dark)

18 – 25°C => 4 hours  
2 – 8°C => 8 hours  
< – 25°C => 14 days

### Acid Phosphatase

(Total and Prostatic)  
18 – 25°C => 4 hours  
2 – 8°C => 8 hours  
< – 25°C => 14 days

- freeze only once

## Warnings and Precautions

Each individual blood donation used for production of SP- MULTICAL. was found to be non-reactive when tested with approved methods for HBsAg, anti-HIV 1+2 and anti-HCV.

As there is no possibility to exclude definitively that products derived from human blood transmit infectious agents, it is recommended to handle the calibrator with the same precautions as valid for patient specimens.

## Preparation

The lyophilisate is vacuum sealed; therefore the vial should be reconstituted carefully with exactly 3 ml of distilled water. Close the vial carefully and allow the calibrator to stand for 30 minutes swirling occasionally. Avoid foaming! Do not shake!

## Procedure

Refer to the package inserts of the reagent kits

### SYMBOLS IN PRODUCT LABELLING

	Authorised Representative		Use by/Expiration Date
	For in-vitro diagnostic use		CAUTION. Consult instructions for use
	Batch Code/Lot number		Manufactured by
	Catalogue Number		(Xi) - Irritant
	Consult instructions for use		Temperature Limitation

## Expected Values

1. The expected value of specific assays are provided on the assay value sheet accompanying each kit, and are lot specific.
2. The expected values are obtained using replicate assay of each manufactured lot of SP-Multical.
3. The individual laboratory values should fall within the expected values.
4. It must however be noted that each laboratory should establish its own normal values and reference range according to GLP.

## Literature

1. Dati F. Reference materials and guidelines for standardization of methods in laboratory medicine. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 1404-26.
2. Moss DW, Henderson AR. Enzymes . In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 2nd ed. Philadelphia: W.B Saunders Company; 1994. p. 735-896.
3. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395)



**Egyptian Company for Biotechnology (S.A.E)**

Obour city industrial area. block 20008 piece 19 A. Cairo. Egypt.

Tel: +202 4489 2248 - Fax: +202 4489 2247

www.spectrum-diagnostics.com

E-mail: info@spectrum-diagnostics.com



**MDSS GmbH**  
Schiffgraben 41  
30175 Hannover, Germany

