

## URI-TRAK

### Test strips for rapid determination of Glucose, Ketones, Protein and pH-value in Urine

REF: 1100 010 (Uri-trak 1) for Glucose  
 REF: 1100 020 (Uri-trak 2) for Glucose & Ketones  
 REF: 1100 030 (Uri-trak 3) for Glucose, Protein & pH  
 REF: 1100 040 (Uri-trak 3 **GK**) for Glucose, Protein & Ketones

#### Intended Use

Spectrum **Uri-trak** is a screening test strips for detection of glucose, ketones, protein and pH-value in urine. Certain configuration of strips may be read instrumentally, using the appropriate Urine Chemistry Analyzers.

#### Principle

**Glucose:** The detection is based on the glucose oxidase-peroxidase-chromogen reaction. Apart from glucose, no other compound in urine is known to give a positive reaction.

**Ketones:** The test is based on the principle of Legal's test. Acetoacetic acid and acetone form with sodium nitroprusside in alkaline medium a violet coloured complex.

**Protein:** The test is based on the "protein error" principle of indicators. The test zone is buffered to a constant pH value and changes colour from yellow to greenish blue in the presence of albumin. Other proteins are indicated with less sensitivity.

**pH:** The test paper contains indicators which clearly change colour between pH 5 and pH 9 (from orange to green to turquoise).

#### Reactive Ingredients

**Glucose:**  
 2% w/w glucose oxidase; 1% w/w peroxidase; 10% w/w TMB; 70% w/w Buffer; 17% w/w nonreactive ingredients.

**Ketones:**  
 5% w/w sodium nitroprusside; 95% w/w Buffer.

**Protein:**  
 0.2% w/w tetrabromophenol blue; 97.4% w/w Buffer; 2.4% w/w nonreactive ingredients.

**pH:**  
 0.2% w/w methyl red; 2.8% w/w bromothymol blue; 97% w/w nonreactive ingredients.

#### Specimen Collection

Use a fresh urine sample that is less than 2 hours old and place it in a clean, dry container. **Do not centrifuge.** The presence of usual urine preservatives will not affect the test results.

#### Instructions for use

Dip the test strip for approximately 1 second into the fresh urine. Draw it across the rim of the container to remove excess urine. After 30 to 60 seconds compare the test strip with the colour scale. The best time for comparison is after 30 seconds. Colour changes that take place after more than 2 minutes are of no significance. When tested the urine should not be older than 2 hours.

SYMBOLS IN PRODUCT LABELLING			
	Authorised Representative		Use by/Expiration Date
	For in-vitro diagnostic use		CAUTION. Consult instructions for use
	Batch Code/Lot number		Manufactured by
	Catalogue Number		(Xi) - Irritant
	Consult instructions for use		Temperature Limitation

#### Evaluation – Sources of Error

**Glucose:** Pathological glucose concentrations are indicated by a colour change from green to bluish green. Yellow or greenish test fields should be considered negative or normal. The colour fields correspond to the following ranges of glucose concentrations: neg. (yellow), neg. or normal (greenish), 50, 150, 500 and 1000 mg/dl or neg. (yellow), neg. or normal (greenish), 2.8, 8.3, 27.8 and 55.5 mmol/l. The influence of ascorbic acid (vitamin C) has been largely eliminated. An inhibitory effect is produced by gentisic acid. Falsely positive reactions can be produced by a residue of peroxide containing cleansing agents.

**Ketones:** The test is more sensitive to acetoacetic acid than to acetone. Values of 10 mg/dl acetoacetic acid or 50 mg/dl acetone are indicated. The colour fields correspond to the following acetoacetic acid values: 0 (negative), 25(+), 100(++) and 300(+++) mg/dl or 0 (negative), 2.5(+), 10(++) and 30(+++) mmol/l. Phenylketones in higher concentrations interfere with the test, and will produce variable colours.  $\beta$ -Hydroxybutyric acid is not detected. Phthalein compounds interfere by producing a red colouration.

**Protein:** The minimum sensitivity of the test strip is 10 mg protein/dl urine. The colour fields correspond to the following ranges of albumin concentrations: negative, 30, 100 and 500 mg/dl or negative, 0.3, 1.0 and 5.0 g/l. Falsely positive results are possible in alkaline urine samples (pH L 9), after infusions with polyvinylpyrrolidone (blood substitute), after intake of medicaments containing quinine and also by disinfectant residues in the urine sampling vessel. The protein colouration may be masked by the presence of medical dyes (e.g. methylene blue) or beetroot pigments.

**pH:** The pH value of fresh urine of healthy people varies between pH 5 and pH 6. The colour scale gives a clear distinction of pH value between pH 5 and pH 9.

#### Reagent Storage and Stability

Store at room temperature (15 to 30 °C) out of direct sunlight. Do not use after expiration date.

#### Precautions

The Kit contain a non-poisonous and harmless desiccant. In case this desiccant is swallowed accidentally, then drink plenty of water.

#### Disposal

Please dispose all used dipsticks in accordance with your local laws and regulations.

ORDERING INFORMATION	
CATALOG NO.	QUANTITY
1100 010 1100 020 1100 030 1100 040	100 Test Strip for each