

Mueller Hinton Agar

REF.	Pack size
1406 001	100 gm
1406 002	500 gm

Intended Use

Mueller Hinton Agar is used for determination of susceptibility of microorganisms to antimicrobial agents isolated from clinical samples

Background

Bauer, Kirby, Sherris and Tuck recommended Mueller Hinton Agar for performing antibiotic susceptibility tests using a single disk of high concentration. Mueller Hinton Agar is mainly used for the primary isolation of Neisseria species. It is specified in FDA Bacteriological Analytical Manual for food testing. The major use of Mueller-Hinton Agar is for Antimicrobial Susceptibility Testing (AST). It has become the standard medium for the Bauer-Kirby method, and is specified by the Clinical & Laboratory Standards Institute (CLSI) formerly the National Committee for Clinical Laboratory Standards (NCCLS) and The European Committee on Antimicrobial Susceptibility Testing (EUCAST).

Principle

Casein acid hydrolysate and beef infusion supply amino acids and other nitrogenous substances, minerals, vitamins, carbon and other nutrients to support the growth of microorganisms. Starch acts as a protective colloid against toxic substances that may be present in the medium. Hydrolysis of starch during autoclaving provides a small amount of dextrose, which is a source of energy.











Components	gm/Liter
Casein Acid Hydrolysate	17.5
Beef Extract Powder	2.0
Starch	1.5
Agar	17.0
Final pH (at 25°C)	7.3 ± 0.1

Preparation, Storage and Stability

Store the dehydrated medium at 10-30°C and use before the expiry date on the label. Store the prepared medium at 2-8°C. After the desired amount of medium is taken out, replace the cap tightly to protect from hydration.

Procedure

1. Suspend 38 g of the powder in 1 L distilled water and mix well.
2. Heat with frequent agitation to dissolve the powder completely.
3. Sterilize by autoclaving at 121°C for 15 minutes.
4. Cool to 45-50°C then mix well and pour into sterile petri plates.

SYMBOLS IN PRODUCT LABELLING		
	Authorised Representative	 Use by/Expiration Date
	For in-vitro diagnostic use	 CAUTION. Consult instructions for use
	Batch Code/Lot number	 Manufactured by
	Catalogue Number	 (Xi) - Irritant
	Consult instructions for use	 Temperature Limitation

Quality Control

Appearance

1-Dehydrated Appearance : beige, homogeneous and free flowing powder.

2-Prepared Appearance : Prepared medium is hazy and light to medium yellow.

3-Cultural Response : Cultural characteristics after 18-24 hours at 35-37°C.

Organisms	Growth
<i>Escherichia coli</i>	Good to luxuriant
<i>Staphylococcus aureus</i>	luxuriant
<i>Enterococci faecalis</i>	luxuriant
<i>Pseudomonas aeruginosa</i>	luxuriant

Interpretation of the results

Refer to appropriate documents for correct zone sizes.

Precautions

Numerous factors can affect results: inoculum size, rate of growth, medium formulation and pH. Strict adherence to protocol is required to ensure reliable results.

Bibliography

1. Mueller and Hinton, 1941, Proc. Soc. Exp. Bio. And Med; 48:330.
2. Bauer et al, 1966, Am. J. Clin. Patho., 45:493.
3. US Food and Drug Adm; 1998, Bacteriological Analytical Manual, 8th Ed; Rev. A, AOAC, International, Gaithersburg, Md.
4. National Committee for Clinical Laboratory Standards. 2000. Approved Standard: M2-A7.
5. Bauer A. W., Kirby W. M., Sherris J. C. and Turck M. (1966) Amer. J. Clin. Path. 45. 493-496.
6. Ryan K. J., Schoenkecht F. D. and Kirby W. M. (1970) Hosp. Pract. 5. 91-100.
7. CLSI Performance Standards for Antimicrobial Disk Susceptibility Testing; CLSI Supplement M100S Wayne PA: Clinical and Laboratory Standards Institute
8. CLSI Performance Standards for Antimicrobial Susceptibility Tests; Approved Standard. CLSI Document M02. Wayne PA: Clinical and Laboratory Standards Institute
9. Pollock H. M. et al. (1986) J. Clin. Microbiol. 24. 1-6.
10. European Committee on Antimicrobial Susceptibility Testing. EUCAST disc diffusion Test Methodology www.eucast.org



Egyptian Co for Biotechnology - Spectrum Diagnostics (S.A.E)
 Obour city industrial area. block 20008 piece 19 A. Cairo. Egypt.
 Tel: +202 4489 2248 - Fax: +202 4489 2247
www.spectrum-diagnostics.com
 E-mail: info@spectrum-diagnostics.com



MDSS GmbH
 Schiffgraben 41
 30175 Hannover, Germany

